

RZ Medizintechnik GmbH  
Unter Hasslen 20 | D-78532 Tuttlingen  
Phone: +49 (0)7462/9470-0  
Fax: +49 (0)7462/9470-50  
sales@rz-medizintechnik.com  
www.rz-medizintechnik.com



RZ  
Medizintechnik  
GmbH  
*Progress in Surgery*

## ENT Shaver System



news  
02.2015

# ENT SHAVER SYSTEM

This shaver system is equipped with the most up to date and powerful micro motor which delivers high levels of reliability. The angled and the straight Shaver Handpiece works as a transmission unit, that transforms the rotating speed of the attached electronic motor in the rate of 4:1.

Advantages:

- Concave and convex shaver blades to prevent "dead" spots
- Integrated cooling system
- 6 shaver blade positions for all working situations

Set consists of:

- Electronic Console 700-100-500
- Foot Pedal 700-100-511
- Micro Motor 700-100-512
- Support for irrigation bag 700-100-515



## Technical Data

Dimensions	290 x 300 x 110 mm
Weight	5.2 kg
Voltage	100 - 240 V AC
Input Frequency	50 - 60 Hz
Pump Capacity	100 ml of fluid per min. at 100% flow rate

252-650-005

ENT Shaver System,  
complete set



252-650-001  
Shaver Handle, angled



252-650-002  
Shaver Handle, straight



Suction control without steps



700-100-515  
Support for irrigation bag  
included



252-650-003  
Tubing Set, sterile, 3m, single use,  
10 pcs. per package



700-100-512  
Micro Motor, 120mm length,  
Ø 20mm, 3m cable length



Irrigation pump unit







700-100-511  
Footswitch, 2m cable length



252-650-004  
Spray for care and maintenance  
of shaver handle


# SHAVER BLADES




Cleaning brush included



Integrated cooling system



Safe and easy locking of each blade position



Shaver blade can be positioned in 6 different positions



Length: 105mm  
5 pcs. per package, non-sterile

**FULL RADIUS RESECTOR**

Order Number	Position	Angle	∅
252-651-025	concave	straight	2.5mm
252-651-026	convex	straight	2.5mm
252-651-035	concave	straight	3.5mm
252-651-036	convex	straight	3.5mm
252-651-040	concave	straight	4.0mm
252-651-041	convex	straight	4.0mm
252-651-135	concave	15°	3.5mm
252-651-136	convex	15°	3.5mm
252-651-140	concave	15°	4.0mm
252-651-141	convex	15°	4.0mm
252-651-235	concave	40°	3.5mm
252-651-236	convex	40°	3.5mm
252-651-240	concave	40°	4.0mm
252-651-241	convex	40°	4.0mm
252-651-335	concave	60°	3.5mm
252-651-336	convex	60°	3.5mm
252-651-340	concave	60°	4.0mm
252-651-341	convex	60°	4.0mm



Length: 105mm  
5 pcs. per package, non-sterile

**AGRESSIVE FULL RADIUS**

Order Number	Position	Angle	∅
252-652-025	concave	straight	2.5mm
252-652-026	convex	straight	2.5mm
252-652-035	concave	straight	3.5mm
252-652-036	convex	straight	3.5mm
252-652-040	concave	straight	4.0mm
252-652-041	convex	straight	4.0mm
252-652-135	concave	15°	3.5mm
252-652-136	convex	15°	3.5mm
252-652-140	concave	15°	4.0mm
252-652-141	convex	15°	4.0mm
252-652-235	concave	40°	3.5mm
252-652-236	convex	40°	3.5mm
252-652-240	concave	40°	4.0mm
252-652-241	convex	40°	4.0mm
252-652-335	concave	60°	3.5mm
252-652-336	convex	60°	3.5mm
252-652-340	concave	60°	4.0mm
252-652-341	convex	60°	4.0mm

# SHAVER BLADES



Length: 105mm  
5 pcs. per package, non-sterile

## TRI SPIKE CUTTER


Order Number	Position	Angle	Ø
252-653-025	concave	straight	2.5mm
252-653-026	convex	straight	2.5mm
252-653-035	concave	straight	3.5mm
252-653-036	convex	straight	3.5mm
252-653-040	concave	straight	4.0mm
252-653-041	convex	straight	4.0mm
252-653-135	concave	15°	3.5mm
252-653-136	convex	15°	3.5mm
252-653-140	concave	15°	4.0mm
252-653-141	convex	15°	4.0mm
252-653-235	concave	40°	3.5mm
252-653-236	convex	40°	3.5mm
252-653-240	concave	40°	4.0mm
252-653-241	convex	40°	4.0mm
252-653-335	concave	60°	3.5mm
252-653-336	convex	60°	3.5mm
252-653-340	concave	60°	4.0mm
252-653-341	convex	60°	4.0mm



Length: 105mm  
5 pcs. per package, non-sterile

## TRIMMER

Order Number	Position	Angle	Ø
252-655-025	concave	straight	2.5mm
252-655-035	concave	straight	3.5mm
252-655-040	concave	straight	4.0mm



Length: 130mm  
5 pcs. per package, non-sterile

## HOODED ROUND BURR

Order Number	Position	Angle	Ø
252-654-055	concave	straight	5.5mm